

Date: 07/21/2016

SEARCH NWL Results reported by Year (N) - Eligible Participants

Enter Name of

Methodology used for each Assay below

(add additional rows as needed)

Enter an 'X' into each box for results that correspond with the assay method

Assay	2002	2003	2004	2005	2006	2007	2008	2009	2010	2011	2012	2013	2014	2015	2016
GAD65	430	1713	1535	1158	935	---	---	---	---	---	---	---	---	---	---
Index	X	X	X	X	X										
GAD65DK	---	---	---	---	---	---	---	6528	553	287	---	1628	173	3209	---
DK								X	X	X	X	X	X	X	
IAA	430	1713	1535	1158	852										
Index	X	X	X	X	X										
IA2	430	1714	1535	1158	935										
Index	X	X	X	X	X										
IA2DK	---	---	---	---	---	---	---	6528	553	287	---	1628	173	3209	---
DK								X	X	X	X	X	X	X	
ZnT8	---	---	---	---	---	---	---	---	---	---	---	---	2158	3256	---
Index												X	X	X	
Tchol	419	1705	1523	1157	967	929	881	877	509	24	1217	1189	718	136	---
Roche 917, Enzymatic	X	X	X	X	X										
Roche Modular P, Enzymatic					X	X	X	X	X	X	X	X	X	X	X
HDLChol	419	1705	1523	1157	967	929	881	877	509	24	1217	1189	717	136	---
Roche 917, Dextran Sulfate, 50 MW precipitation	X	X	X	X	X										
Roche Modular P, Dextran Sulfate, 50 MW precipitation					X	X	X	X	X	X	X	X	X	X	X
LDLChol	419	1705	1523	1157	967	929	881	877	509	24	1217	1189	717	136	---
TG <400mg/dL Derived (Friedewald Equation), TG ≥ 400mg/dL Full Beta Quantification	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
VLDLChol	419	1705	1523	1157	967	929	881	877	509	24	1217	1189	717	136	---
TG <400mg/dL Derived (Friedewald Equation), TG ≥ 400mg/dL Full Beta Quantification	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Trig	419	1705	1523	1157	967	929	881	877	509	24	1217	1189	718	136	---
Roche 917, Enzymatic	X	X	X	X	X										
Roche Modular P, Enzymatic					X	X	X	X	X	X	X	X	X	X	X
HbA1c	415	1699	1547	1170	1001	1028	980	949	550	25	1233	1214	732	143	---
BioRad Variant II	X	X	X	X	X	X									
Tosoh G-7						X	X	X							
Tosoh G-8 Start Date: 1/16/09								X	X	X	X	X	X	X	X
Albumin (Urine)	413	1630	1453	1129	947	933	882	823	477	19	1048	1017	644	126	---
Siemens reagents, BNII Nephelometer	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Creatinine (Urine)	413	1630	1453	1129	947	933	882	823	477	19	1048	1017	644	126	---
Roche 917, Jaffe Method	X	X	X	X	X										
Roche Modular P, Jaffe Method					X	X	X								
Roche Modular P, Enzymatic Start Date: 8/12/08							X	X	X	X	X	X	X	X	X
Albumin/Creatinine Ratio (Urine)	374	1445	1270	976	851	839	783	745	426	18	923	899	569	116	---
Calculation	X	X	X	X	X	X	X	X	X	X	X	X	X	X	X
Glucose	417	1710	1523	1157	968	929	881	877	510	24	1217	1189	718	136	---
Roche 917, Hexokinase	X	X	X	X	X										
Roche Modular P, Hexokinase					X	X	X	X	X	X	X	X	X	X	X
Glucose Tolerance Test															
Glucose0Min	6	115	454	600	168	18	5	5	9	---	---	---	---	---	---
Glucose30Min	5	113	449	590	306	176	81	60	57	---	---	---	---	---	---
Glucose60Min	6	113	447	592	310	174	79	59	57	---	---	---	---	---	---
Glucose90Min	6	113	18	---	---	---	---	---	---	---	---	---	---	---	---
Roche 917, Hexokinase	X	X	X	X	X										
Roche Modular P, Hexokinase					X	X	X	X	X						

